

biopen® | Q2

Auto MTS Electroporation

- EMS
- MTS
- LED



Q2 is a beauty instrument that integrates microneedling, EMS electroporation, and LED skin care. It is suitable for beauty salons and home-use. Equipped with an induction base, it stops when the pen is placed and starts when it is off the base. It can meet the needs of skin care, hair growth and beard growth etc.

Focus on microneedling therapy in beauty and cosmetics industry


Facial introduction, auxiliary absorption

- Own Professional R&D team
- 100,000-class dust-free workshop
- Brand owner and original manufacturer of well-known brand Dr.pen



QUALIFICATION CERTIFICATE

- ISO 13485 Medical devices quality management standard
- EU CE/RoHS Certification
- US FDA Certification
- US FCC Certification
- US EMC Certification
- EU LVD Certification
- Bio Pen US/EU Trademark Certificate
- Dr. Pen US/EU Trademark Certificate



Selling point

Triple effects Rejuvenate the skin

Microneedling
Accelerate skin repair

EMS Electroporation
Promote skin collagen synthesis

LED skin care
Calm and soothes skin



25 Pole Pulse EMS Electroporation*

Stable and average energy
Three-dimensional cell stimulation

EMS electroporation micro needle
Stimulate cells and accelerate collagen synthesis



EMS Microcurrent energy field
The ingenious arrangement of positive and negative poles allows each adjacent needle to form countless combinations of network forms, thus forming a dense microcurrent energy field.

* 25 pole pulsed EMS electroporation is 25 independent needles arranged ingeniously to complete a current quickly and transmit the current evenly.
* The needle can be replaced with a 9-pin needle, in this case it is a 9-pin pulse electroporation, and so on.

Four skincare in one use

Microneedling+EMS Electroporation+Blue Light+Red Light

Repair skin by Microneedling
Microneedles create thousands of micro channels to initiate skin repair.


EMS Electroporation
EMS microcurrent enters skin through the channel created by microneedling and achieves a multiplier effect.

470nm Blue Light
Reduce acne, cleanse, calm and soothe
Repair sensitive skin and balance oil

630nm Red Light
Promote collagen regeneration, shrink pores
Increase skin elasticity, brighten skin tone



For Skin Care



- ✓ Fine lines, wrinkles
- ✓ Pore size, acne scars
- ✓ Skin discoloration
- ✓ Helps improve skin texture
- ✓ Stretch mark on belly
- ✓ Large areas of acne
- ✓ Stretch mark on thighs


For Hair Growth & Beard Growth

It's not that the hair restorer has no effect, but your deliver method is not ideal.

Say goodbye to traditional skin care methods.



- Create the micro channels and deliver hair growth liquid directly to the deep layer in order to increase the absorption efficiency.
- EMS stimulates cells to produce ATP, providing more energy for hair growth while strengthening hair roots.
- Increase absorption.
- Strengthen hair and prevent hair loss.
- Shorten the resting period of the hair follicle.



Advantages

Pause when it is placed and starts when off

Comes with an induction base. The pen immediately pauses work when placed on the base, while remaining the current speed level, and continues to work with the same speed when off the base.

Induction Base Pause, Place and Display

Low battery warning

Q2 has low battery warning to help users know the remaining battery power intuitively. It can better control the working hours.

Smart touch screen

One touch button, Clear and convenient.

Product Parameters

Specification			
Brand	bio pen	Model	Q2
Power Model	Wireless	Handle Material	ABS
Speed Level	5 Speeds	Available Cartridges	9 Pins 25 Pins
EMS microcurrent	3 levels	Battery capacity	700 mAh
EMS frequency	0.5-0.7 Hz	Adapter	US/EU/UK/AU/KR
LED light	Blue light, Red light	G.W.	194*146*67 cm/1 pc/ 676 g
Needle Adjustment Range	0-1.5mm	1 carton	41*52*30 cm / 20 pcs / 15 kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box ② bio pen Q2 ③ 25 pins * 2 pcs ④ Manual ⑤ Inductive base ⑥ Adapter ⑦ Cable